

MINIMUM HARDWARE REQUIREMENTS



OS: MACOS 13 Ventura,
14 Sonoma, 15 Sequoia

PROCESSOR: M2

DISPLAY: 1920x1080 (Full HD)

INTERNET: stable DSL connection

MacBook PRO 16"
iMac Series



OS: Windows 10, Windows 11

PROCESSOR: Intel i5, Ryzen 5

RAM: minimum 16 GB

HARD DISK: minimum 300 GB

GRAPHIC CARD: 4GB or higher vRAM

To ensure better AI performance, use the NVIDIA
GTX or QUADRO line (compute capability \geq 5.0)

DISPLAY: 1920 x 1080 for optimal visualization

INTERNET: stable DSL connection



BROWSER: Google Chrome

INTERNET: stable DSL connection

*In certain jurisdictions, including the United States and the European Union, RealGUIDE is a registered trademark of 3Diemme S.r.l.. The use of the ® symbol is limited to the jurisdictions where the trademark is registered.